

Patient monitoring solutions for your medical facility

Monitors • Anesthesia • Electrosurgical • Surgery Lights • Surgery Tables and much more



Waveline EZ
Portable Patient Monitor



Waveline Touch
Patient Monitor



Waveline Pro
O.R. Monitor



Envoy TS
Central Monitoring Station



Transport | Emergency



Bedside | Vital Signs



O.R. | Critical Care



Central Station Monitoring

When monitoring multiple parameters & multiple patients, we have you covered.

Envoy TS

Monitor patients throughout your facility from a single location

The configurable Envoy TS Central Monitoring Station includes wireless ECG transmitters. No other system provides the high level of flexibility offered by the Envoy TS. Each channel is independent, allowing the use of stationary, wireless transport and ECG telemetry transmitters on a single system.

Parameters: ECG monitoring using portable transmitters, or get complete integration with DRE-brand monitors for ECG, respiration, temperature, SpO₂ and NIBP.



NEW

On-site software upgrades available on DRE Waveline Series Monitors



Waveline EZ

Portable Patient Monitor with an intuitive touch-screen

At only 6 pounds, the Waveline EZ is one of the most mobile multi-parameter monitors on the market. It features an 8" high-resolution touch-screen that provides immediate operation and a VGA connector for use with LCD monitors.

Parameters: ECG, resp., SpO₂, NIBP and temp. Available with EtCO₂ and a printer.

Designed for: Triage, emergency environments, ambulance use, and patient transport.

Waveline Touch

Lightweight Patient Monitor with a large touch-screen

The Waveline Touch features a 12.1" color TFT display with a touch-screen interface, multi-lead ECG monitoring, ST segment monitoring and 72 hours trend data analysis.

Parameters: ECG, resp., SpO₂, NIBP, heart rate and temp. Available with dual IBP, EtCO₂, built-in recorder and a VGA connector for use with LCD monitors. Also available with Optional Plug-In Anesthesia Gas Module.

Designed for: Pre-operation and post-operation, ER, transport, CCU, and conscious sedation rooms.

Waveline Pro

O.R. Monitor available with EtCO₂ and anesthesia gas monitoring

The lightweight Waveline Pro displays eight colorful waveforms on a crisp 15" high-resolution screen. It works with EMR systems, and the optional agent bench is part of the main unit—making it a very compact anesthesia monitor.

Parameters: ECG, resp., dual temp., SpO₂ and NIBP. Features Available with EtCO₂, Dual IBP, anesthetic gas monitoring, cardiac output monitoring and a printer.

Designed for: General surgery, ICU, and anesthetic agent monitoring.

Envoy Telemetry Transmitter

Portable ECG transmitter for Envoy TS Telemetry System

Flexible, durable Envoy TS telemetry transmitters feature user-selectable lead wire configurations and provide an extended battery life.