

Full Functionality, Intensity and Versatility



The Xavier-C10 Portable LED Headlight provides physicians with a high intensity 10-watt LED mobile light source that is designed specifically for the operating room. The comfortable, lightweight design and long-life battery pack are good for hours of continuous use.

Physicians wanting to supplement their existing surgical light set-up in the operating room will find that the Xavier-C10 is one of the most powerful, comfortable and cost-effective LED headlight systems on the market today.

Features:

- High Intensity 10-watt LED bulb rated for 25,000 - 50,000 hours of use
- 8,100 Foot Candles at distance of 30cm
- 5,000 Foot Candles at distance of 50cm
- 6,000° Kelvin color temperature
- Clear field of illumination prevents hot-spots and shadows
- Rechargeable battery pack guarantees 4.5 hours of maximum intensity light
- Charge on the go with built-in belt clips on each battery pack
- LED Battery Life Meter to measure remaining time before recharging
- Easily-accessible On/Off switch and intensity controls
- Pinpoint spot adjustment
 - » 7cm - 1.25cm spot at distance of 30cm
 - » 9cm - 2cm spot at distance of 40cm
- Lightweight, adjustable headband for user comfort
- Sterilizable positioning bar adjusts from right to left, up and down

Standard Configuration:

- 10-watt portable LED headlight
- Battery pack — 2 rechargeable packs with each system
- 100 - 240V Battery charger
- Convenient carrying case

