

## European designed hydraulic delivery bed with safe, simple and convenient operation.

### Features:

- Simple and convenient manual hydraulic operation
- Constructed of quality 304 stainless steel
- Easy to maintain and sterilize

### Standard Accessories include:

- Anesthesia Screen
- Shoulder Support
- Body Support
- Arm Support
- Leg Support
- Foot Plate



### Technical Data:

Tabletop (Length/Width)	69.29" L x 23.62" W
Tabletop Max. Elevation	35.83" H
Tabletop Min. Elevation	25.98" H
Trendelenburg/ Reverse Trendelenburg	12°/22°
Back Plate Adjustment	+45°/-10°