

Power and long life on the go, in a battery operated suction unit.

Performance

No need to worry about power when you need it. This high performance unit delivers a vacuum range of 2-21 inHg. Caregivers can easily adjust the vacuum with the aspirator's standard, easy-read gauge.

Long Lasting

The high performance battery allows 1.5 hours of run time, twice that of comparable units. A "Low Battery Light", standard in this unit, informs the caregiver that the unit needs to be re-charged thus preventing interruptions in use.

Convenient

No need to buy a special canister. The DM-530 Portable Aspirator uses any standard 800 cc canister and is visible in the Carry Bag while unit is in operation. It is re-chargeable in its convenient, specially designed Carry Bag, which also stores the charger, Suction Tubing and optional Auto DC Power Cord.

Durable

The ABS plastic cabinet is durable, attractive, easy-to-clean and infection control friendly. There are no open areas on the unit to harbor bacteria and dirt.

Protected

A dual pump protection system features two safeguards. A float shut-off in the canister to protect overflow and an Inline Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pump.

Dependable

The DM-530 Portable Aspirator carries a full 2 year warranty on the unit and a 6 month warranty on the battery.



DM-530

Portable Aspirator

Equipment for the way *you* operate

Specifications:

Dimensions (in Carry Bag):	12" (L) x 9 1/2" (W) x 10 3/16" (H)
Weight (with Carry Bag and all Accessories):	9.6 lbs. (4.4 kg)
Shipping Weight:	11 lbs. (5 kg)
Electrical:	
Aspirator unit:	12 Volts DC Battery
Charger/Power Supply:	100-240 Volts AC
Operating Temperature Range:	50°F to 104°F (10°C to 40°C)
Storage Conditions:	
Temperature:	-4°F to 104°F (-20°C to 40°C)
Humidity:	Max 95% non-condensing
Vacuum Range:	50 mmHg to 535 mmHg (2 inHg to 21 inHg)
Vacuum Buildup:	0 to 300 mmHg in 3.1 sec. (with 800cc Suction Canister and AC Power Cord connected) 0 to 300 mmHg in 3.5 sec. (with 800cc Suction Canister connected and running on Battery)
Suction Canister:	800cc Disposable with Float Shut-off
Power Consumption:	20 watts
Flow:	14 lpm (unrestricted flow, without suction canister, running on AC power)

