

E-ST Series

Ultrasonic Pocket Doppler

Equipment for the way *you* operate

Easy to use, compact ultrasonic Doppler



Compact Li-ion Battery Charger



Large Backlit LCD Display



Built-in Audio Recorder



Interchangeable Waterproof Probe



Additional features include:

- Accurate FHR detection with clear sound
- Supports 2MHz, 3MHz, 4MHz and 8MHz waterproof probes
- Automatic and manual FHR counting
- Up to four hours of continuous use
- Automatic power-off

Specifications

LCD display	45 mm x 25 mm
Heart rate range	50~210 bpm
Resolution	1 bpm
Audio sampling frequency	4KHz/8KHz
Recording length	480s/240s
Audio output power	0.5W
Auto-shutoff	1 minute
Size	32 mm (D) x 85 mm (W) x 138 mm (H)
Weight	290±5g (includes one battery)

Model	FHR Display	Built-in Audio Recorder	Backlit Display	Probe Auto Detection Display	Battery Status Display	AC Recharge
E-ST Lite	No	No	No	No	No	No
E-ST Basic	Yes	No	No	Yes	No	No
E-ST Basic A	Yes	No	Yes	Yes	Yes	No
E-ST Pro	Yes	Yes	Yes	Yes	Yes	No
E-ST II	Yes	No	Yes	Yes	Yes	Yes
E-ST Pro II	Yes	Yes	Yes	Yes	Yes	Yes
Options	Rechargeable Li-ion battery, charger, earphone, carry bag, 9V alkaline battery, adapter, probe (2MHz, 3MHz, 4MHz, 8MHz)					