

The ideal anesthesia machine for today's high-traffic operating room

The Integra SL3 is expertly designed with clinician-focused features including customizable flowmeter display options, a standard mounting system and a fully integrated ventilator and absorber. The Selectatec-compatible backbar can fit two to three vaporizers.



< Vital signs monitor sold separately

Features:

- Choose which flow control visualization display is right for you:
 - » Standard/cascade flow tubes
 - » Cascade flow tubes with electronic flow display: numeric view of gas flow
 - » Virtual Flow Display (VFD): numeric and graphic view of flow tubes with virtual spinning bobbins and traditional total flow tube
- Selectatec-compatible backbar holds up to three vaporizers
- Hospital Information System (HIS) capabilities when using electronic flow display or VFD.
- Convenient pull-out writing space for note-taking and charting
- Fully-integrated CO2 absorber and bellows unit with ventilator interface
- Standard GCX™ compatible mounts
- Large capacity drawers for easy storage



Standard Configuration:



AV-S Ventilator: A multi-functional model that includes three advanced spontaneous support modes and electronic PEEP. The AV-S Ventilator is designed for both adult and pediatric use, and includes integrated oxygen monitoring and spirometry, plus seamless integration with the anesthesia machine and A200SP absorber.



A200SP Absorber: Includes ventilator interface that provides mode switching triggered by the bag or ventilator control. The A200SP features a quick-release canister for easy removal and cleaning. All parts are autoclavable.



Sigma Delta Vaporizer: Available with multi-choice fill system options.

Specifications:

- Overall Dimensions: 28" (71cm) W x 27.5" (70cm) D x 55" (139cm) H
- Work Surface Height: 34" (86cm)
- Work Surface: 15" (38cm) x 10" (25cm) | holds 66 lbs. (30 kg) — evenly distributed
- Writing Tablet: 12" (30cm) x 8.5" (22cm) | 22 lbs. (10 kg) — evenly distributed
- Top Shelf: 28" (71cm) x 13.5" (35cm) | 66 lbs. (30 kg) — evenly distributed
- Drawers: 4.5" (12cm) x 21.5" (54.5cm) x 13.5" (35cm) | 22 lbs. (10 kg) — evenly distributed
- Castors: 5" (125mm) diameter — all with brakes
- Weight: 275 lbs. (125 kg)

