

## Versatile surgical table that provides stability and the positioning you need — including horizontal sliding

### Features:

- Electro-hydraulic operation and a sleek design.
- Surface is made of 304 stainless steel.
- Safe and reliable with simple, convenient controls.
- Automatic base lock by remote control.
- One button repositioning — great for X-ray and C-arm use.
- Seamless integration into orthopedic environments.



### Specifications:

- 770 lbs. / 350 kg. weight capacity
- 550 lbs. / 250 kg. articulated weight capacity
- Tabletop Length/Width: 1,970 mm, 500 mm (77.5 in, 19.7 in)
- Tabletop Elevation: High: 880 mm, Low: 600 mm (34.6 in, 23.6 in)
- Trendelenburg/ Reverse Trendelenburg: 30°
- Lateral tilt: 18°
- Head plate adjustment: Up: 90°, Down: 60°
- Leg plate adjustment: Up: 15°, Down: 90°, Outward: 90°
- Back plate adjustment: Up: 90°, Down: 45°
- Kidney bridge: 120 mm (4.7 in)
- Horizontal slide: 350 mm (13.8 in)



Smooth patient positioning by remote control



Standard Accessories:

