


Universal power procedure table

Designed for diagnostics, taking blood-samples, gynecology, dermatology and minor surgery.

- Fully electric multi-function table with four actuators and standard foot pedal control.
- The electronic control box and the actuators are versatile and easy to service. Featuring an automatic reset switch that returns the table to its starting position.
- Available in 110 and 220 V cc. configurations. The galvanized steel frame is fully covered with ABS for an attractive appearance as well as safety.
- The standard configuration includes an anatomic mattress, adjustable and removable headrest, two guides for medical accessories and a towel roll-holder.
- The anatomic mattress is made of top-quality foam that includes a five-year warranty. The upholstery is designed for medical use, long-lasting and easy to clean. Plus it is resistant to water, air and disinfectants.

Optimized for

- General Medical
- Dermatology
- Podiatry
- Plastic, Cosmetic and Aesthetics
- OB/GYN
- ENT/Otolaryngology
- Oral & Maxillofacial
- Dental
- Orthodontics



Shown with optional medical arm-rests.





Shown with optional accessories.

Color Options



Optional Accessories

1. I.V. Rod
2. Curtain-holder
3. Medical Arm-rest
4. Jointed Blood Sample Taking Arm-rest
5. Jointed Thigh Support
6. Hand-set Control
7. Fixed Clamp
8. Rotating Clamp
9. Colposcope Support



Table Dimensions

