


## Maximum Efficiency and Safety



The new **PREMIO** electric hospital bed is the perfect solution for hospitals, surgery centers and teaching facilities. This innovative hospital bed is fully adjustable for maximum ease of use and comfortable procedure. With control systems for both the patient and the caregiver, this bed provides optimal safety and prevention of disease.

### **Features:**

- » Safe and firm construction from steel with powder coating and platform of ABS plastic with ventilation holes
- » Head and foot board easily removable
- » The drop down side rails each have a patient control panel for Hi-Lo, back, knee and auto-chair adjustments with an additional exterior control panel for caregivers to use for hygienic operation
- » Foot-end attendant control panel for easy control of patient position
- » Auto-regression function
- » Locking control panel prevents unexpected bed movements
- » Mattress Stopper
- » I.V. Pole Holes and Drainage Bag Hangers: Six I.V. pole holders are located around the bed and hangers for drainage bags on each side of the bed
- » Auto-Chair Function simultaneously manipulates the back and knee rests for comfortable positioning
- » Back-up battery



**Specifications:**

**Safe working load:** 496 lbs (225 kg)

**Castor Diameter:** 5"

**Motor Standard:** IP 54 (optional IP 66)

**Equipment Standard:** IEC60601-2-38

**Back Rest CPR Quick Release Lever:** optional

**Central Break and Steering:** optional

**Dimension:**

- » Overall Bed Length: 85.5in (2170mm)
- » Overall Bed Width: with side rails up- 42in (1070mm), with swivel side rails up- 44in (1112mm), with side rails stowed- 38 in (960mm)

» Overall Mattress Base Length: 76in (1940mm)

» Overall Mattress Base Width: 36.5in (925mm)

» Overall Bed Height: 18-29.5 in (450-750 mm)

**Angle:**

» Back rest: 0-75 deg

» Knee gatch: 0-22 deg

**Trendelenburg & Reverse Trendelenburg:** Optional

**Back rest angle indicator:** Yes

**Auto-chair:** Yes

