

The **Torino EXL** automatic operating table is an excellent choice for use in a variety of operations. A micro touch remote control operates the quiet hydraulic lift system, flexible adjustments of the head, back, and seat plates, and the built-in kidney bridge. The **Torino EXL** features long horizontal sliding, C-arm and x-ray access.



#### Features

- C-arm and X-ray interface capabilities to meet imaging requirements
- Micro touch remote control with low noise movements
- Flexible adjustments on head plate, back plate, seat plat and built-in kidney bridge
- Quiet and reliable hydraulic lift system
- Stainless steel design for excellent stability



### Technical Specifications

- Length/width: 2070 mm/550 mm; 81.5 in/21.7 in
- Elevation (down/up): 700 mm/1,000 mm; 27.6 in/39.4 in
- Trendelenburg: 25°
- Lateral tilt (left, right): 15°
- Head plate (up/down): 45°/90°
- Leg plate (up/down/outward): 15°/90°/90°
- Back plate (up/down): 75°/20°
- Kidney bridge up: 110 mm; 4.3 in
- Horizontal sliding: 300 mm; 11.8 in
- Power: 220 V, 50 Hz, 1 Kw or 110 V, 60 Hz, 1.5 Kw

### Standard Accessories

- Anesthesia screen: 1 piece
- Shoulder support: 1 pair
- Body support: 1 pair
- Arm support: 1 pair
- Leg support: 1 pair
- Foot plate: 1 pair

