

True ECG 12T

Twelve-channel ECG

Equipment for the way *you* operate

Features

- 8.4" color touch-screen
- Alphanumeric keyboard and one touch operation
- Sleep mode to save energy and extend LCD life
- Automatic measurement and interpretation tested with authoritative CSE database
- Real-time waveform freezing
- 120 seconds ECG waveform reviewable
- Two-step exercise test with periodic recording
- Built-in rechargeable Li-ion battery
- 200 ECGs in internal memory



True ECG 12T

Available with
PC-based
Data Management Software/
Smart ECG Viewer



- Manual Measurement
- Rhythm Viewing
- Waveforms Contrast
- Auto Measurement Result
- Multiple Report Formats

Excellent cost performance with built-in high-resolution thermal printer

- Preview report for both real-time printing and file storage
- Average template waveform with accurate interpretation tested with CSE DATABASE, which supplies measurement, diagnosis, and Minnesota code
- Doctor can edit the printing style, printing speed, sensitivity, and used filter to get better waveforms instead of repeating patient examination
- Freeze and review waveform timely to grasp every arrhythmia
- Internal storage including waveforms, measurement and diagnosis which could also be enlarged by flash disk or PC hard disk.



3111 N. University Drive, Suite 625
Coral Springs, Florida 33065.
Telephone: 954-757-1250
Email: info@accessiblemediequip.com

www.accessiblemediequip.com

Accessible MediEquip (USA)