

The Versailles P100 mobile surgery table goes beyond the O.R.

TRANSPORT



SURGERY



The Versailles P100 is a complete surgery solution that can **move patients from the pre-operative area to the operating room and then to the recovery area** — with no need to transfer the patient from a stretcher to a surgery table and back.

RECOVERY



The Versailles P100 does it all, increasing throughput by incorporating **stretcher mobility** into a **multi-purpose surgery table!**

PRE-OP



The complete multi-purpose mobile surgery table

Powered patient trolley/operating table is a total solution for hospitals and surgery centers.

The Versailles P100 is more than a surgery table.

It is a complete surgery solution that can move patients from the pre-operative area to the operating room and then to the recovery area — with no need to transfer the patient from a stretcher to a surgery table and back. The Versailles P100 does it all by incorporating the mobility of a stretcher into a multi-purpose surgery table.

Features

- » Four remote-activated powered functions
- » Lateral tilt
- » Excellent stability
- » Full C-arm coverage
- » Superior surgical access
- » Traversing patient platform
- » 550-lb. (250-kg) weight limit



Pictured with optional armboards



Accessories

Surgical Headrests

Armboards

Leg and Foot Positioners

plus much more



Full C-arm coverage

Stretcher maneuvering is incorporated into the complete surgery system

Stretcher Mobility

Patient safety rails fitted for patient comfort and safety.

Retractable steering mechanism engaged for maximum maneuverability.



Gynecology features

Patient platform traversed to foot end for access.

Head section removed for anesthesia access.

Leg section removed for access.

Fitted with optional powered leg supports.



Shoulder Arthroscopy features

Patient platform traversed to foot end for access.

Optional arthroscopy backrest fitted.



Ophthalmology

Patient platform traversed to head end for access.

Ophthalmic head section fitted.



± 12° Lateral Tilt



Removable head section



Push button rotating cot-sides



Push button release and self-locking leg section

Technical Specifications

Height (min)	23.8" (605 mm)
Height (max)	39.5" (1,005 mm)
Trolley length	82.7" (2,100 mm)
Side Rail Width	25.8" (655 mm)
Patient Safety Rail Width	30.3" (770 mm)
Brake Width (Brakes off)	33" (835 mm)
Mattress Length	80.3" (2,040 mm)
Mattress Width	23.6" (600 mm)
Mattress Depth	3" (75 mm)
Base Length	47" (1,190 mm)
Trendelenburg	20°
Reverse Trendelenburg	12°
Lateral Tilt	+12°/-12°
Backrest Articulation	0 - 65°
Head Section Articulation	+30°/-30°
Head End Traverse	9" (230 mm)
Foot End Traverse	9" (230 mm)
C-arm Access	100%
Max. Patient Weight Limits	550 lbs (250 kg)
Trolley Weight	352 lbs (160 kg)
Lightweight Leg Section Weight	13.2 lbs (6 kg)

Specifications subject to change without notice.

Trolley Features

- » Powered Functions
- » Lateral Tilt
- » Head End Traverse
- » Foot End Traverse
- » K8 Pressure Care Mattress
- » Removable Ultra-Lightweight Leg Section
- » Removable and Rotating Cot-sides
- » Push Button Release/Self Locking Controls
- » 5th Wheel Steering Mechanism
- » Advanced Four-Caster Braking System
- » E Size Oxygen Cylinder Trough (Side Loading)
- » Transfusion Pole with Mounting Position
- » 150 mm Castors

Optional Features

- » Spare Transfusion Pole
- » Sculpted Mattress for Head Support
- » Ophthalmic Head Support
- » Removable Shoulder Arthroscopy Backrest
- » Padded Patient Safety Rail Covers - Gray
- » Padded Patient Safety Rail Covers - Child Print
- » Articulating Removable Leg Section
- » Removable Leg Section with Full Length Side Rail
- » Press Mattress Complete
- » **Accessory Trolley: Complete storage solution**
- » Accessory Trolley Spare Trough
- » Spare Powered Trolley Hand Control
- » Spare Powered Trolley Charging Cable
- » CD Side Cylinder Support Bracket



Removable ultra lightweight leg section



K8 pressure care mattress

