

Versailles M75

Mobile Surgery Table (Manual)

Equipment for the way *you* operate

The complete multi-purpose **mobile operating table**

Manual mobile surgery table is a total solution for hospitals and surgery centers.



- Mechanically powered functions
- Increased surgical access
- Lightweight and maneuverable
- 75% C-arm access
- Low cost solution
- 350-lb. (160 kg) weight limit

Transport and Recovery

- Cotsides fitted for patient comfort and safety
- Unique retractable steering mechanism engaged for maximum maneuverability

Gynecology features

- Lightweight leg section removed for access
- Head section removed for anesthesia access
- Fitted with optional powered leg supports



ENT/Dental

- Patient positioned on standard articulating head section



Low trolley height

- Helps patients more easily get on and off the trolley



Removable head section



Push button rotating cotsides



Push button release and self leg section



Removable ultra lightweight leg section



K8 pressure care mattress

Versailles M75

Mobile Surgery Table (Manual)

Equipment for the way *you* operate

Technical Specifications – Manual System

Height (min)	22.6" (575 mm)
Height (max)	36.8" (935 mm)
Trolley length	82.7" (2,100 mm)
Side Rail Width	25.8" (655 mm)
Cotside Width	30.3" (770 mm)
Brake Width (Brakes off)	33" (835 mm)
Mattress Length	80.3" (2,040 mm)
Mattress Width	23.6" (600 mm)
Mattress Depth	3" (75 mm)
Base Length	41.5" (1,055 mm)
Trendelenburg	20°
Backrest Articulation	0 - 80°
Head Section Articulation	+30°/-30°
C-arm Access	75%
Max. Patient Weight Limits	352 lbs (160 kg)
Trolley Weight	275 lbs (125 kg)
Lightweight Leg Section Weight	13.2 lbs (6 kg)

Trolley Features

- K8 Pressure Care Mattress
- Removable Ultra Lightweight Leg Section
- Removable and Rotating Cotsides
- Push Button Release and Self Locking Controls
- 5th Wheel Steering Mechanism
- Advanced 4 Castor Braking System
- E Size Oxygen Cylinder Trough (Side Loading)
- Transfusion Pole with Mounting Position
- 150 mm Castors

Optional Features

- Spare Transfusion Pole
- Sculpted Mattress for Head Support
- Ophthalmic Head Support
- Removable Shoulder Arthroscopy Backrest
- Padded Cotside Covers - Gray
- Padded Cotside Covers - Child Print
- Articulating Removable Leg Section
- Removable Leg Section w/ Full Length Side Rail
- Press Mattress Complete
- Accessory Trolley
- Accessory Trolley Spare Trough
- Spare Powered Trolley Hand Control
- Spare Powered Trolley Charging Cable
- CD Side Cylinder Support Bracket



Optional Accessory Trolley

- Complete storage solution for Versailles M75 accessories and other surgery items



C-arm access

- Provides 75% C-arm access

Manual positioning maximizes affordability

