

Low energy, low heat, high efficiency

Using about 50% less energy than similar halogen lighting, the Vista LED exam/diagnostic light is compact, maneuverable, and provides great light output with a large flood pattern. By incorporating LED technology, the Vista LED features the best light output, color temperature and heat output (extremely low) in its class. The Vista LED is cost effective and represents a leap in technology and efficiency.

Features:

- Large, flood pattern — approximately 6" (150 mm) diameter.
- 50,000 lux of homogenous light output at 1 meter.
- 4,300° Kelvin color temperature at 24" (61 cm) distance.
- Life rating of 50,000 hours.
- Uses about 50% of the energy of similar traditional halogen lighting.
- Extremely cool output — significantly less heat generated than halogen lighting.
- IR absorbing, color-correcting glass filter.
- Available in wall, portable floor, single or dual ceiling mounts, and combination ceiling mounts with Vision or Vortex lights.



Superior Articulation:

28" (711 mm) articulating arm with integral compression spring.

- 23" extension arm.
- Arm rotates 40° up and down from horizontal plane.
- Lamphead and arm both rotate 540°.
- Positioning or optional sterilizable handle mounted on perimeter of lamphead.

Choose from a variety of configurations

