

## Customizable illumination for your operating needs

The Vista Pro 3 LED Examination Light is a full-featured, customizable light ideal for use in ICU, dermatology and neonatal environments. Adjust between 3 distinct color temperatures and 5-stage dimming for the perfect illumination in every setting.

### Features:

- Brilliant 40,000 Lux output at 0.5 meters
- 5-stage dimming
- Alternate between 3 preset color temperatures: 4,400° Kelvin, 3,800° Kelvin, 3,000° Kelvin
- Long-life LED bulbs use less energy than halogen alternatives
- Available in spring arm or gooseneck configuration
- Light arm rotates 360° at the base

### Specifications:

- Light head diameter: 4.7" (12cm)
- Weight: 3.1 lbs (1.4 kg)
- Light body: polycarbonate
- Power: 18 W
- Voltage: 100-240V
- Frequency: 50/60Hz



Power and dimming control panel located on the light head



### Available in multiple options:

Rail, table or wall mount with spring arm or gooseneck configurations or portable stand.

